

# MEMORY TASK APPLICATION - INSTRUCTIONS

What's Included:

You should have the following files/folders inside the extracted folder:

- main.exe                <- The application you will run
- images/                <- Contains all experimental images
- result\_output/        <- Where participant data will be saved
- README.txt            <- This instruction file

How to Run the Application:

1. Double-click 'main.exe' to launch the task.
2. A dialog box will appear asking for:
  - Subject ID (e.g., sub01)
  - Session Number (e.g., 1)
  - Task Mode (choose "self\_paced" or "fixed")
3. If a data file for that subject/session already exists, the program will ask whether to overwrite it.
4. The memory task will then begin. The participant will see a series of images and respond using the keyboard:
  - Press B if the item is NEW
  - Press N if the item is OLD
  - Press M if the item is SIMILAR
5. After completing the task, the program will automatically return to the initial screen to allow entry for the next participant.

## Where to Find the Experimental Data:

All session data are saved in the 'result\_output/' folder under each participant's subfolder.

For example:

result\_output/

|-- sub01/

    |-- session\_1\_self\_paced.csv

## CSV File Format Explanation:

Each '.csv' file contains a row for each trial and includes the following columns:

- trial : Trial number
- image : Image filename shown
- type : Trial type ('old', 'similar', or 'new')
- expected : Expected correct response label ('N', 'O', or 'M')
- response : Actual key pressed by participant ('N', 'O', 'M', or 'F' for no response)
- correct : Whether the response was correct (1 = correct, 0 = incorrect, 'F' = no response)

Tip: Always double-check that a '.csv' file was created after each session!

## Questions or bugs?

Contact the developer if you have any issues or questions running the application.